

Installation Completion Confirmation

Dentium
For Dentists By Dentists

Product Name	<input type="checkbox"/> rainbow TM CT : <input checked="" type="checkbox"/> 2 in 1 <input type="checkbox"/> 3 in 1 <input type="checkbox"/> bright TM CT : 1T (<input type="checkbox"/> 2 in 1 <input type="checkbox"/> 3 in 1) 2T (<input type="checkbox"/> 2 in 1 <input type="checkbox"/> 3 in 1) 2T ST (<input type="checkbox"/> 2 in 1 <input type="checkbox"/> 3 in 1) <input type="checkbox"/> rainbow TM MCT		
Office Name	ARIDENTIST CLINIC	Office Name	
Address	ROVINARI	Address	ROMANIA
Serial Number	H22 AXD0124	Serial Number	
Measure Tool	Digital Multimeter (0224)		

No.	Inspection	Status and Operation	Result
1	Exterior	1 Check if there is any abnormality in the appearance of the equipment - After installation, check the appearance of the equipment for damage.	<input checked="" type="checkbox"/> OK <input type="checkbox"/> NG
2	Power	1 Check input voltage (110~240V) - The voltage of the wall power source to which the equipment power is applied.	[240.6] V
3	Operation	1 Check if the LED changes from red to green when power is applied	<input checked="" type="checkbox"/> OK <input type="checkbox"/> NG
		2 Check normal operation of Column	<input checked="" type="checkbox"/> OK <input type="checkbox"/> NG
		3 Check the laser operation	<input checked="" type="checkbox"/> OK <input type="checkbox"/> NG
		4 Check the normal operation of the equipment origin control	<input checked="" type="checkbox"/> OK <input type="checkbox"/> NG
		5 Check the operation of temple support	<input checked="" type="checkbox"/> OK <input type="checkbox"/> NG
4	Workstation	1 Check the PC running status	<input checked="" type="checkbox"/> OK <input type="checkbox"/> NG
		2 Check the network connection status	<input checked="" type="checkbox"/> OK <input type="checkbox"/> NG
5	Software	1 Check the installation ① rainbow ImageViewer ② 3D Viewer ③ rainbow Ceph (3 in 1 only)	<input checked="" type="checkbox"/> OK <input type="checkbox"/> NG
6	Performance	1 rainbow ImageViewer File, DB connection status check	<input checked="" type="checkbox"/> OK <input type="checkbox"/> NG
		2 Check the operation status of the exposure switch LED in the exposure ready state	<input checked="" type="checkbox"/> OK <input type="checkbox"/> NG
		3 Check the status indicator while exposure ① exposure switch: Green → Red ② status light : Green → Yellow	<input checked="" type="checkbox"/> OK <input type="checkbox"/> NG
		4 image output and storage	<input checked="" type="checkbox"/> OK <input type="checkbox"/> NG
7	Accessory	1 Check that you have received the components that match the specifications you ordered: ① Standard/edentulous chinrest ② edentulous bite ③ sinus chinrest ④ Tooth Bite ⑤ Temple support ⑥ Product manual ⑦ Carpus jig (3 in 1 only)	<input checked="" type="checkbox"/> OK <input type="checkbox"/> NG
8	User training	1 Check that you have been trained on how to use the product	<input checked="" type="checkbox"/> OK <input type="checkbox"/> NG

Installer :	IONUT RUBOIAN,	(Signature)
Sales representative :	IULIAN FLOREA	(Signature)
Customer :	ARIDENTIST CLINIC	(Signature)

Customer Satisfaction Survey

Dentium
For Dentists By Dentists

This survey will be used as the basic information to improve user convenience, product quality, and source of new development.

Customer	AR DENTIST CLINIC	Country	ROMANIA
Address	ROVINARI	Phone	0040728083186
Product	Rainbow CT	Date	11.09.2025

1. What was the most important consideration when selecting our product? (Multiple)

- | | |
|--|---|
| <input checked="" type="radio"/> 1) Price | <input checked="" type="radio"/> 2) Image quality |
| 3) Brand | 4) Product Quality |
| <input checked="" type="radio"/> 5) Customer service | 6) Others () |

2. Select the most satisfying answers using our product. (Multiple)

- | | |
|---|---|
| <input checked="" type="radio"/> 1) Image quality () | 2) Image viewer () |
| 3) Convenience () | <input checked="" type="radio"/> 4) Service () |
| 5) Others () | |

3. Please answer each of the following questions. (If 'Poor or lower', describe your opinions below)

	Excellent	Good	Average	Poor	Very Poor
1) Product design	<input checked="" type="radio"/> 1	<input checked="" type="radio"/> 2	3	4	5
2) Convenient image acquisition	<input checked="" type="radio"/> 1	2	3	4	5
3) Convenient software manipulation	<input checked="" type="radio"/> 1	2	3	4	5
4) Satisfied with the Pan image quality	<input checked="" type="radio"/> 1	2	3	4	5
5) Satisfied with the CBCT image quality	<input checked="" type="radio"/> 1	2	3	4	5
6) Satisfied with the Ceph image quality	<input checked="" type="radio"/> 1	2	3	4	5
7) Easily understand the user's manual	<input checked="" type="radio"/> 1	2	3	4	5

Others



4. Please tell us your comments or complaints about our products or service in general.

