

Installation Completion Confirmation

Dentium
For Dentists By Dentists

| | | | |
|---------------|--|---------------|--|
| Product Name | <input type="checkbox"/> rainbow TM CT : <input type="checkbox"/> 2 in 1 <input checked="" type="checkbox"/> 3 in 1 <input type="checkbox"/> bright TM CT : 1T (<input type="checkbox"/> 2 in 1 <input type="checkbox"/> 3 in 1) 2T (<input type="checkbox"/> 2 in 1 <input type="checkbox"/> 3 in 1) 2T ST (<input type="checkbox"/> 2 in 1 <input type="checkbox"/> 3 in 1) <input type="checkbox"/> rainbow TM MCT | | |
| Office Name | SC DENTAL HOPE SRL | Office Name | |
| Address | Strada Henry Barbusse 59-61, Cluj-Napoca | Address | |
| Serial Number | L02DXD0624 | Serial Number | |
| Measure Tool | Digital Multimeter (25542) | | |

| No. | Inspection | Status and Operation | Result |
|-----|---------------|--|--|
| 1 | Exterior | 1 Check if there is any abnormality in the appearance of the equipment - After installation, check the appearance of the equipment for damage. | <input checked="" type="checkbox"/> OK <input type="checkbox"/> NG |
| 2 | Power | 1 Check input voltage (110~240V) - The voltage of the wall power source to which the equipment power is applied. | [230] V |
| 3 | Operation | 1 Check if the LED changes from red to green when power is applied | <input checked="" type="checkbox"/> OK <input type="checkbox"/> NG |
| | | 2 Check normal operation of Column | <input checked="" type="checkbox"/> OK <input type="checkbox"/> NG |
| | | 3 Check the laser operation | <input checked="" type="checkbox"/> OK <input type="checkbox"/> NG |
| | | 4 Check the normal operation of the equipment origin control | <input checked="" type="checkbox"/> OK <input type="checkbox"/> NG |
| | | 5 Check the operation of temple support | <input checked="" type="checkbox"/> OK <input type="checkbox"/> NG |
| 4 | Workstation | 1 Check the PC running status | <input checked="" type="checkbox"/> OK <input type="checkbox"/> NG |
| | | 2 Check the network connection status | <input checked="" type="checkbox"/> OK <input type="checkbox"/> NG |
| 5 | Software | 1 Check the installation ① rainbow ImageViewer ② 3D Viewer ③ rainbow Ceph (3 in 1 only) | <input checked="" type="checkbox"/> OK <input type="checkbox"/> NG |
| 6 | Performance | 1 rainbow ImageViewer File, DB connection status check | <input checked="" type="checkbox"/> OK <input type="checkbox"/> NG |
| | | 2 Check the operation status of the exposure switch LED in the exposure ready state | <input checked="" type="checkbox"/> OK <input type="checkbox"/> NG |
| | | 3 Check the status indicator while exposure ① exposure switch: Green -> Red ② status light : Green -> Yellow | <input checked="" type="checkbox"/> OK <input type="checkbox"/> NG |
| | | 4 image output and storage | <input checked="" type="checkbox"/> OK <input type="checkbox"/> NG |
| 7 | Accessory | 1 Check that you have received the components that match the specifications you ordered. ① Standard/edentulous chinrest ② edentulous bite ③ sinus chinrest ④ Tooth Bite ⑤ Temple support ⑥ Product manual ⑦ Carpus jig (3 in 1 only) | <input checked="" type="checkbox"/> OK <input type="checkbox"/> NG |
| 8 | User training | 1 Check that you have been trained on how to use the product | <input checked="" type="checkbox"/> OK <input type="checkbox"/> NG |

| | | |
|------------------------|--------------------|-------------|
| Installer : | IULIAN FLOREA | (Signature) |
| Sales representative : | FLORIN SCHEOPU | (Signature) |
| Customer: | SC DENTAL HOPE SRL | (Signature) |

